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Pre-configured Small Variant Demo Filter (29)

Pre-configured TSO 500 RNA Demo Filter (0)


Pre-configured Focal CNV Demo Filter (1)

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


	Benign			Oncogenic			
Evidence Strength	Very Strong	Strong	Supporting	Supporting	Moderate	Strong	Very Strong
POINTS	-8	-4	-1	+1	+2	+4	+8
Population Data	MAF is >5%	MAF is >1%		Absent in population databases			
Functional Data		Well-established functional studies show no oncogenic effects				Well-established functional studies supportive of an oncogenic effect	
Predictive Data			Silent mutation (no predicted impact on splicing)		Missense change at an amino acid residue where a different missense change determined to be oncogenic has been documented	Same amino acid change as a previously established oncogenic mutation	Null variant in tumor suppressor
Cancer Hotspots				Cancer hotspots with low frequency of recurrence	Cancer hotspot with moderate frequency of recurrence	Cancer hotspot with high frequency of recurrence	
Computational Evidence			All utilized lines of computational evidence suggest no impact of a variant	All utilized lines of computational evidence support oncogenicity			

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ESCAT

ESMO Scale for Clinical Actionability of Molecular Targets

Summary

Patient ID: 529893
Tissue ID: T25-17898
Report date: 1-Feb-2024

Patient and sample details

Hospital	University Medical Centre
Patient number	529893 (F)
Date of birth	15-May-1967
Tissue ID	T25-17898
Date of sampling	12-Jan-2024
Sampling anatomical site	Liver
Sample type	Biopsy (FFPE)
Histological cancer type	Adenocarcinoma (NSCLC) (PA-report: T24-12811)
Tumour cell content	Pathology-based: 60% Sequencing-based: 53%
Sample ID(s)	M25-02999, M25-03001
Matched normal DNA	M23-0157
Source of normal	Blood

History / Reason for testing

Patient (female) with metastatic NSCLC and a known *ALK* fusion (*EML4-ALK*). Treated with Crizotinib. Sample obtained upon progression after a short initial response. Resistance mechanisms, other targets?

Summary of most relevant findings

Lung adenocarcinoma metastatic biopsy showing:

- *ALK* (*EML4* ex13 – *ALK* ex20) gene fusion. This fusion has also been detected in a previous tumour sample of this patient. No additional (resistance) mutation are observed. Possible indication for second/third-line ALK inhibitor (e.g. Alectinib, Brigatinib, Lorlatinib, Ceritinib)
 - *BRCA1* (p.Ala1708Glu, p.Gln1756fs) inactivation, ESCAT III in lung cancer (possibility for drug repurposing of PARP inhibitors in a clinical trial setting)
 - Homologous recombination deficiency (HRD) (score 0.95) confirming the complete (bi-allelic) inactivation of *BRCA1*.
- One of the observed *BRCA1* (p.Gln1756fs) mutations is detected in the patient's germline DNA and previously reported as Ashkenazi founder pathogenic variant. Referral to a genetics specialist can be considered.
- The combination of two *BRCA1* mutations (bi-allelic) and high a HRD score indicates functional loss of homologous recombination repair function.
- To translate these findings into treatment recommendations, further interpretation of the results should be performed by a physician considering the patient's full medical history, if appropriate with support of a molecular tumour board.*
- An overview of all detected cancer-associated somatic genomic events and germline variants can be found in the report.*

The data on which this report is based has been generated using tests that are performed under ISO-17025:2017 TESTING accreditation. Assay details, sensitivity/specificity specs and clinical validation reports are available via genomics-facility.org/ISO-val-reports/. Based on a minimal tumour purity of at least 20%, the test has a sensitivity of >95% for detection of genomics variants, >95% for copy number alteration and >95% for detection of gene fusions.

Genes/regions that are covered by the assay(s) design and/or that are included in the clinical reporting can be downloaded here: www.genomics-facility.org/assays/. Further details about the bioinformatics pipeline and analysis parameters can be found on GitHub: (<https://github.com/seqfacility/>)

This report has been generated on 1-Feb-2025 using the pipeline version 2.1. A user manual to assist in interpretation of this report can be downloaded here: www.genomicsfacility.org/test-manuals/.

For questions, feedback or complaints please contact support@genomics-facility.org

Results

Patient ID: 529893
Tissue ID: T25-17898
Report date: 1-Feb-2024

Mutations

Assay/analysis QC: X

Gene	Type	Transcript	Position	Variant	Effect	Reads	VAF	LoH	Germline	Hotspot	Class
<i>TP53</i>	TSG	NM_000546.6	17:7577124	c.814G>A (p.Val272Met)	missense	372 / 610	61%	yes	no	yes	Likely pathogenic
<i>BRCA1</i>	TSG	NM_007294.4	17:41209079	c.5266dup (p.Gln1756fs)	frameshift	297 / 571	53%	no	yes	yes	Pathogenic
<i>BRCA1</i>	TSG	NM_007294.4	17:41215920	c.5123C>A (p.Ala1708Glu)	missense	231 / 662	35%	no	no	yes	Pathogenic
<i>TERT</i>	ONCO TSG	NM_198253.3	5:1295228	c.-124C>T	promoter	238 / 366	65%	no	no	yes	Pathogenic
<i>KEAP1</i>	TSG	NM_203500.2	19:10602583	c.996del (p.Gly332fs)	frameshift	287 / 521	55%	yes	no	no	Likely pathogenic
<i>KDM6A</i>	ONCO TSG	NM_001291415.2	X:44922943	c.1804G>A (p.Gly602Arg)	missense	199 / 398	51%	no	no	no	VUS

Copy number alterations

Assay/analysis QC: X

Gene	Type	Transcript	Region	Event	Gene copy number	Chr arm copies
<i>RB1</i>	TSG	NM_000321.3	13q14.2	loss	0	4
<i>NCOA2</i>	ONCO	NM_006540.4	8q13.3	gain	10	2

Gene fusions

Assay/analysis QC: X
Comment Warning, low alignment score*

Gene 5'	Transcript 5'	Exon 5'	Gene 3'	Transcript 3'	Exon 3'	Phasing	Reads
<i>EML4</i>	NM_019063.5	13	<i>ALK</i>	NM_004304.5	20	In-frame	135/211 (88%)
<i>HEY1</i>	NM_012258.4	3	<i>NCOA2</i>	NM_006540.4	10	exon-skipping	6/14 (43%)

*Gene fusions show a relatively low number of aligned reads, indicating a lower sensitivity/specificity, likely due to lower quality of the input tissue. Caution is required for interpretation of the results

Mutational signatures

Assay/analysis QC: X

Signature	Score	Threshold [range]	Result
Tumour mutation burden (TMB)	2.4	10 [0-100+]	Low
Microsatellite (in)stability	0.1	4 [0-100+]	MSS (stable)
Homologous recombination deficiency	0.95	0.5 [0-1]	Deficient

Genomics aberrations are reported based on the HG37 genome assembly and the MANE select transcript IDs (or the transcript indicated) (last update 01-Jan-2023). For reported (>3) variants, a loss-of-heterozygosity status has been determined to indicate complete bi-allelic loss of the wildtype allele. Variants are classified as (likely) pathogenic, variant of unknown significance (VUS). (Likely) benign variants are not included in this report. Variants that are (likely) present in the patient's germline are indicated. Genes are marked as Oncogenes (ONCO) and/or tumour suppressor genes (TSG) and variants in known hotspots codons are indicated. (source: Jackson CKB, as per 1-Jan-2024). More details and an overview of the genes that are considered for variant reporting can be found here: www.genomics-facility.org/assays/variant-reporting/

Copy number alterations are reported for copy gains x chromosomal arm copies). For copy losses, only events are reported that show a complete loss (0 copies). Gene losses are reported for TSG only, gene gains (amplifications) are called for ONCO. More details and an overview of the genes that are considered for fusion reporting can be found here: www.genomics-facility.org/assays/CNV-reporting/

Gene fusion detection has been performed using RNA-based sequencing. Fused gene products are reported in case of in-frame events, including in-frame events that require exon skipping. Gene fusions are preferentially detected on the MANE select transcripts but can also include alternative transcripts. The number and percentage of fusion reads are reported. More details and an overview of the genes that are considered for fusion reporting can be found here: www.genomics-facility.org/assays/fusion-reporting/

Reported genomic signatures include the tumour mutational burden (TMB) as number of missense muts/Mb, microsatellite instability (MSI), and homologues recombination deficiency (HRD). Details of the used procedure and bioinformatics pipeline can be found on GitHub (<https://github.com/seqfacility/>)

Biomarkers clinical evidence

Patient ID: 529893
Tissue ID: T25-17898
Report date: 1-Feb-2024

Clinical actionability annotation

Gene/Biomarker	Alteration	Match level	Treatment type/ drug	ESCAT	Source
<i>ALK</i>	<i>ALK-EML4 fusion</i>	fusion	Alectinib, Lorlatinib, Crizotinib, Ceritinib, Brigatinib	I	CKB-JAX
<i>ALK</i>	<i>ALK-EML4 fusion</i>	fusion	Alectinib, Lorlatinib, Crizotinib, Ceritinib, Brigatinib	I	OncokB
<i>BRCA1</i>	p.Ala1708Glu, p.Gln1756fs	inactivation	Olaparib, Talazoparib, Niraparib	III	CKB-JAX
<i>BRCA1</i>	p.Ala1708Glu, p.Gln1756fs	inactivation	Rucaparib	III	CKB-JAX
<i>BRCA1</i>	p.Ala1708Glu, p.Gln1756fs	oncogenic	Olaparib, Niraparib, Rucaparib	III	OncokB
HRD	Positive	signature	PARP inhibitors	III	CKB-JAX

Clinical trial biomarker matching

Gene/Biomarker	Alteration	Match level	Trial name/number	Phase	Source
<i>BRCA1</i>	p.Gln1756fs	inactivation	NCT04826341 NCT03845296 NCT05700721	I/II II	CKB-JAX
HRD	Positive	signature	NCT04826341 NCT05700721	I/II II	CKB-JAX

Findings for potential follow-up

Gene/Biomarker	Alteration	Follow-up	Remark
<i>BRCA1</i>	c.5266dup (p.Gln1756fs)	Referral to genetic counselling/cancer genetics specialist	Pathogenic variant also observed in the patient germline DNA. Reported as Ashkenazi founder pathogenic variant.
<i>HEY1::NCOA1</i>	<i>HEY1</i> ex3 – <i>NCOA1</i> ex10	Confirmatory testing is required if considered for the pathological diagnosis.	Low RNA reads, potential passenger due to <i>NCOA2</i> amp.

Reported genomic events are matched against knowledgebases (source: Jackson CKB, OncokB, as per 1-Jan-2024) for associated clinical evidence regarding potential treatment approaches. ESCAT levels are used for clinical relevance; only evidence items with an ESCAT tier I, II or III are reported. Clinical evidence matching is performed for (likely) pathogenic variants (VUS are not included), copy numbers alterations, fusions and genomics signatures. Only the evidence items with the highest tier match are shown.

Clinical trial matching is performed based on the genomics event and tumour type information only. Other clinical trial eligibility criteria are not considered, and need to be evaluated further if inclusion in a clinical trial is considered.

Shrnutí

	SNV/InDels	SVs	Del/Amp	MSI	TMB	HRD
FISH	✗	✓	✓	✗	✗	✗
PCR	✓	✓	✗	✓	✗	✗
NGS	✓	✓	✓	✓	✓	✓

	SNV/InDels	SVs	Del/Amp	MSI	TMB	HRD
PCR (RT-PCR)	✗	✓	✗	✗	✗	✗
NGS	✗	✓	✗	✗	✗	✗

	SNV/InDels	SVs	Del/Amp	MSI	TMB	HRD
IHC	✓*	✓*	✓*	✓	✗	✗

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